

In the Specification

Please substitute the following paragraph on page 26, beginning at line 6:

The beads were allowed to settle for 2 h prior to siphoning off the depleted supernatant. The settled beads were then transferred to 50 ml Falcon tubes (Becton Dickinson Biosciences, UK) and recovered by centrifugation at 200g. The beads were transferred to a 30 ml column (Biorad, UK) and then washed with 300 mls of cold phosphate-buffered saline, 0.5 M NaCl, pH 8.3. The column was eluted with 0.1 M citric acid, pH 3.0 into 10 x 2 ml fractions in glass tubes containing 0.4 ml of 2.75M Tris, pH 8.5 for immediate neutralization of the citric acid. Fractions containing CD28TFc according to analysis by 12% sodium dodecyl sulphate polyacrylamide gel electrophoresis (SDS-PAGE) were concentrated to 0.5ml or no higher concentration than 20 mg/ml using a Centriprep 10 concentrator (Millipore Corp). The protein was then applied to a Superdex 200 H/R gel filtration column (Amersham Biosciences) pre-equilibrated with 10 mM Hepes, 140 mM NaCl, pH 7.4 (HBS buffer) for up to three successive runs. Eluting fractions were monitored for absorbance at 280 nm. Protein-containing fractions were examined by SDS-PAGE. Each cycle of batch purification yielded ~9 mgs of CD28TFc; up to 7 ~~batch-purifications~~purifications was required to deplete all the CD28TFc per set of 6 cell factories (~5 litres of starting tissue culture supernatant).